



Precision Dietary IgA (Saliva) - KIT Collection & Contents

Collection Requirements

Refrain from eating or drinking anything for 20 minutes before collection.

Shipping Requirements

ONLY Ship samples Monday—Thursday

Kit Contents

Contact Us if any items are damaged or missing

Quantity	Item Description
1	Requisition form
1	Biohazard bag with gauze pad
1	Prepaid shipping label
1	Laboratory shipping pack
4	Salivette saliva collection tubes

Sample Stability

Refrigerate samples at 2—8 C° until shipment. We ask that samples are received within 6 days of collection. Do not freeze the samples.

Contact Us

9 Dunwoody Park, Dunwoody, GA 30338
(678) 736-6374 | info@precisionpointdiagnostics.com
www.precisionpointdiagnostics.com



Precision Dietary IgA (Saliva) - KIT Collection & Contents

Collection Requirements

Refrain from eating or drinking anything for 20 minutes before collection.

Shipping Requirements

ONLY Ship samples Monday—Thursday

Kit Contents

Contact Us if any items are damaged or missing

Quantity	Item Description
1	Requisition form
1	Biohazard bag with gauze pad
1	Prepaid shipping label
1	Laboratory shipping pack
4	Salivette saliva collection tubes

Sample Stability

Refrigerate samples at 2—8 C° until shipment. We ask that samples are received within 6 days of collection. Do not freeze the samples.

Contact Us

9 Dunwoody Park, Dunwoody, GA 30338
(678) 736-6374 | info@precisionpointdiagnostics.com
www.precisionpointdiagnostics.com

Collection Instructions

1. Collect the saliva in the each of the tubes provided.
2. Saliva is to be collected in physician's office at time of appointment.
3. Remove the cap from the tube and remove the cotton plug. Do not remove the suspended insert.
4. Place the cotton plug into your mouth and gently chew the cotton for 1-2 minutes until fully saturated. Replace the cotton back into the Salivette.
5. Repeat the process for all Salivettes.
6. Using a permanent marker, write the patients first and last name, date of birth, and date of collection on each sample tube.
7. Store sample in the refrigerator at 2-8 °C until ready to ship.

Shipping Instructions

1. Make sure the specimen tubes are properly labeled.
Improperly labeled tubes cannot be processed.
2. Place the tube(s) inside of the biohazard bag and seal the bag. Do not remove gauze pad.
3. Complete the Test Requisition Form, fold and place the form within the outer pocket of the biohazard bag.
4. Place the biohazard bag into kit box
5. Place the kit box inside the lab mailing pack.
6. Attach prepaid shipping label.
7. Take the package to the shipping carrier's storefront location or call the carrier to schedule a pickup.

Contact Us

9 Dunwoody Park, Dunwoody, GA 30338
(678) 736-6374 | info@precisionpointdiagnostics.com
www.precisionpointdiagnostics.com

Collection Instructions

1. Collect the saliva in the each of the tubes provided.
2. Saliva is to be collected in physician's office at time of appointment.
3. Remove the cap from the tube and remove the cotton plug. Do not remove the suspended insert.
4. Place the cotton plug into your mouth and gently chew the cotton for 1-2 minutes until fully saturated. Replace the cotton back into the Salivette.
5. Repeat the process for all Salivettes.
6. Using a permanent marker, write the patients first and last name, date of birth, and date of collection on each sample tube.
7. Store sample in the refrigerator at 2-8 °C until ready to ship.

Shipping Instructions

1. Make sure the specimen tubes are properly labeled.
Improperly labeled tubes cannot be processed.
2. Place the tube(s) inside of the biohazard bag and seal the bag. Do not remove gauze pad.
3. Complete the Test Requisition Form, fold and place the form within the outer pocket of the biohazard bag.
4. Place the biohazard bag into kit box
5. Place the kit box inside the lab mailing pack.
6. Attach prepaid shipping label.
7. Take the package to the shipping carrier's storefront location or call the carrier to schedule a pickup.

Contact Us

9 Dunwoody Park, Dunwoody, GA 30338
(678) 736-6374 | info@precisionpointdiagnostics.com
www.precisionpointdiagnostics.com